

Short Working Distance Sealed Lens Barrel

- **>0.5 mm working distance probe barrel**

- **For turbid solutions or solid samples**



- **Performs** better than ball probe lens barrel



Part Number	Description
SWD-PTFE-6	Short working distance probe barrel



Comparison of the three Probes when measuring Polystyrene and NaAc



- Results are normalized. NaAcetate signal drops as %polystyrene increases. The standard probe signal drops much faster as compared to either the SWD or Ball probes.
- NaAcetate, with no polystyrene added has a normalized signal of ~30,000. At 0.5% polystyrene the signal drops to 5000 counts for the std probe while the SWD probe has a signal of ~18,000 and the BallProbe ~15,000 counts.
- The PS intensity does not track linearly for any probe but even less for the standard probe while the SWD and Ball probes track in a similar fashion, indicating that the effective working distance of these two probes is similar.

Comparison between standard, short working distance, and ball probes

Where the focus is on the SWD lens barrel