

Immersion Lens Barrel

- **Compression sealed immersion lens barrels**

- **Material** of the lens barrel can be adjusted to account for high temperature applications
- **Temperature Range** goes from -40°C up to 300°C (can go over 300°C with gold seal)
- **Pressure Rating** can withstand up to over 3000 psi for high pressure use



Specifications

Barrel Material	316 stainless steel standard (hastelloy, inconel, titanium and other metals available)
Tube Diameter	1/4" to 1"
Tube Length	3", 6" or up to 12" max
O-Ring Seal	Polymer (kalrez, teflon) or metal (gold)
Pressure Rating	Ambient to >3000 psi
Temperature Range	-40°C to 300°C (Kalrez o-ring)
Part Number	SPI-3/6/Custom